

NeuroPAL ID Keyboard Shortcuts

mouse click = select a neuron

mouse double-click = add or remove a neuron

space bar (when cursor in User ID field) = save user ID

The following shortcuts apply when you have already clicked on the image:

0 = toggle the statistical atlas of neuron locations on/off

1 = toggle the red channel on/off

2 = toggle the green channel on/off

3 = toggle the blue channel on/off

4 = toggle the white (panneuronal) channel on/off

5 = toggle the DIC channel on/off

6 = toggle the GFP (reporter) channel on/off

← = decrease the Z slice (go left or dorsal)

→ = increase the Z slice (go right or ventral)

↑ = zoom in

↓ = zoom out

p = toggle image panning on/off

r = restore the image center with no zoom

h = zoom in and toggle through vertical image halves

j = zoom in and toggle through vertical image thirds

k = zoom in and toggle through vertical image quarters

a = toggle auto-naming neurons

b = toggle showing neuron birth times

i = toggle the user IDs on/off